



Microfluidics™

Superior Knowledge | Superior Results



LM20 Microfluidizer™ High Shear Fluid Processor

LM20 Microfluidizer™ High Shear Fluid Processor


Digitally Controlled, Electrically Powered Lab Unit for Small Sample Processing

Microfluidizer™ technology efficiently converts fluid pressure into shear forces, leading industry performance standards in high pressure processing.

A unique solution to maintaining consistent process pressure ensures 100% of your material gets exactly the same treatment. Whether you are working with small-scale lab batches or production volumes, the Microfluidizer processor is unmatched in submicron particle/droplet size reduction, cell disruption, product yield, and guaranteed process scale up.



High efficiency cell disruption with minimal protein denaturation



Recommended For:

- ◆ Emulsions
- ◆ Dispersions
- ◆ Liposomes
- ◆ Cell Disruption
- ◆ Fine Particle Deagglomeration

Unique Benefits of the LM20

- ◆ Achieve unmatched particle/droplet size reduction or cell disruption performance at lower process pressures
- ◆ Precise temperature control
- ◆ Higher product yields after sterile filtration with tighter particle size distribution
- ◆ Save on development time in pilot/production scale with linear volumetric scale-up
- ◆ Enhanced repeatability with easy-to-use digital pressure control
- ◆ Compact design to fit in a fume hood
- ◆ Ensure dependable process performance over time with maintenance reminder and operator alerts



LM20 Microfluidizer™ High Shear Fluid Processor

Operating Principle

The LM20 Microfluidizer processor contains an intensifier pump designed to supply the desired pressure at a constant rate to the product stream. As the pump travels through its pressure stroke, it drives the product at constant pressure through the precisely defined fixed-geometry microchannel within the Interaction Chamber™.

As a result, the product stream accelerates to high velocities, creating shear rates within the product stream that are orders of magnitude greater than any other conventional means. All of the product experiences identical processing conditions, producing the desired results, including: uniform particle and droplet size reduction (often submicron), deagglomeration and high-yield cell disruption.

A removable cooling coil and a cooling bath are incorporated into the design to promote optimal temperature control.

Standard Features


- ◆ 300 ml glass reservoir
- ◆ Diamond Interaction Chamber assembly (Ceramic Interaction Chamber for LM20-20 model)
- ◆ Digital pressure control
- ◆ Integrated stroke counter with maintenance reminder
- ◆ Multiple languages (selectable by user)

Specifications

Pressure Range	Up to 30,000 psi (2068 bar)*
Minimum Sample Size	30 ml
Flowrate Range	Up To 90 ml/min
Feed Temperature	Maximum of 165°F (75°C)
Electrical Power Requirement	208-230/460V, 60Hz, 3PH (190-220/380-415V, 50Hz, 3PH)
Pump Actuator	Electro-hydraulic
Noise Level	<80dB TWA
Dimensions	23"H x 33"W x 26"L (58cm x 84cm x 66cm)
Weight	250 lbs (114 kg)

Options

- ◆ Auxiliary Processing Module™ (APM™)
- ◆ Cooling coil, bath and recirculation assembly
- ◆ Heat exchanger
- ◆ Process pressure gauge
- ◆ Larger capacity glass or stainless steel reservoirs
- ◆ Solvent-resistant gasket materials
- ◆ Seal quench fittings
- ◆ Nitrogen adapter assembly





Microfluidics International Corporation
90 Glacier Drive, Suite 1000,
Westwood, MA 02090, USA

Tel: 617-969-5452 • 800-370-5452
Fax: 617-965-1213

Email: mixinginfo@idexcorp.com
Web: www.microfluidics-mpt.com

Industry Leading Insights & News:
<https://www.microfluidics-mpt.com/blog>

